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

# COMPRESSOR NEBULIZER ULAIZER HOME CN-02 MY

## USER'S MANUAL

## COMPRESSOR NEBULIZER ULAIZER HOME CN-02 MY INSTRUCTIONS

Index .....	
Introduction .....	18
Important Safeguards .....	19
Product Identification .....	21
Usage .....	22
Cleaning .....	23
Disinfection .....	24
Replacement Parts & Optional Accessories .....	24
Troubleshooting.....	26
Specification.....	27
EMC Information.....	28

### 1. INTRODUCTION

Thank you for purchasing the **Ulaizer HOME CN-02 MY**. It is a compact medical device designed to efficiently deliver physician prescribed medication to the bronchial lung passages. With proper care and use, it will provide you with many years of reliable treatment.

This product is developed for the successful treatment of asthma, allergies and other respiratory disorders. It creates a stream of air that travels through clear tube to the nebulizer. When air enters the nebulizer, it will convert the prescribed medication into aerosol mist for easy inhalation.

Your **Ulaizer HOME CN-02 MY** should be used under the supervision of a licensed physician and/or a respiratory therapist. We encourage you to thoroughly read this guidebook to learn about the features of this product.

Any use of this product other than its intended use should always be avoided.

### 2. IMPORTANT SAFEGUARDS

 **Note:** Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

 **Caution!** Failure to read and observe all precautions could result in personal injury or equipment damage.

#### Product cautions

1. To avoid electrical shock: Keep the unit away from water  
Do not immerse the power cord or the unit in liquid  
Do not use while bathing  
Do not reach for a unit that has fallen into water – immediately unplug the unit.
2. Never operate the unit if it has any damaged parts (including power cord), if it has been dropped or submersed in water. Promptly send it to a service center for examination and repair.
3. The unit should not be used where flammable gas, oxygen or aerosol spray products are being used.
4. Disconnect the unit from the electrical outlet before cleaning, filling and after each use.
5. Do not use other attachments unless recommended by manufacturer.
6. This device can be used on conscious persons only.

#### Operating Cautions

1. Connect this product to an appropriate voltage outlet for your model.
2. Do not run this product unattended.
3. Never operate if this unit has a damaged cord or plug, if it has been dropped, or in any way, into water, if it does not work properly. Return it to a service center for repair.
4. If any abnormality occurs, discontinue use immediately until the unit has been examined and repaired.
5. Always unplug the product immediately after use.
6. Never block the air openings of the main unit or place it where the air openings may be obstructed.

### Storage Cautions

1. Do not store the unit under direct sunlight, high temperature or humidity.
2. Keep the unit out of reach of small children.
3. Always keep the unit unplugged while not in use.

### Cleaning Cautions

1. Do not immerse the unit in water. It may damage the unit.
2. Disconnect the unit from the electrical outlet before cleaning.
3. Clean all necessary parts after each use as instructed in this guidebook.

### Important


The compressor motor has a thermal protector, which will shut off the unit before the unit is overheated. When the thermal protector shuts the unit off, please:

1. Switch off the unit.
2. Unplug the unit from the electrical outlet.
3. Wait 30 minutes for the motor to cool down before another treatment. Make sure the air openings are not obstructed. When treatment is finished, shut off the unit and unplug it from the electrical outlet.

When treatment is finished, shut off the unit and unplug it from the electrical outlet.

### 3. PRODUCT IDENTIFICATION

#### ULAIZER HOME CN-02 MY\*

\*Appearance of compressor slightly different from image on box

#### 4. USAGE

1. Remove your ULAIZER HOME CN-02MY unit from the packaging, place it on a flat and stable surface. Be sure that you can easily reach the controls when you are seated and ensure that the ventilation slits are free.
2. Make sure the switch in the "OFF" position.
3. Gently twist the top part of the nebulizer bottle counter-clockwise to disassemble the nebulizer bottle (Figure 1).
4. Fill the nebulizer bottle with the medication prescribed by your physician (Figure 2).

**!** **Note:** Fill with medication directly into the nebulizer bottle, do not fulfill. The recommended fill level is 5ml.

Use medication only on the advice of a physician, and check the appropriate inhalation duration and quantity for your needs.

5. Gently twist the top part clockwise to re-assemble the nebulizer bottle. Be sure the two sections fit well (Figure 3).
6. Connect the air tube with nebulizer bottle (Figure 4).
7. Connect the other end of the air tube to the tube connector on the nebulizer unit by lightly turning it (Figure 5).
8. Firmly connect the mask, mouth piece, nose piece with nebulizer bottle (Figure 6).
9. Hold the nebulizer bottle as shown on the picture and follow the instructions of your doctor or respiratory therapist. Exhale as normal through the mouth piece, mask or nose (Figure 7, 7a).
10. Plug the power cord into an appropriate electrical outlet. Make sure at this stage, the power switch is at "OFF" position.
11. Press power switch to begin your prescribed treatment.
12. Switch off the nebulizer device after treatment. Unplug the power cable from the electrical outlet.


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6

Fig. 7

Fig. 7a

#### 5. CLEANING

**!** **Important**

Before you clean the unit, always switch it off, unplug it and let it cool off.

The nebulizer bottle and the used accessories such as mouthpiece, mask, etc. must be washed with hot water after each use.

Carefully dry all components using a soft cloth. Once all components are fully dry, reassemble the nebulizer bottle and place the components in a dry, sealed container or carry out disinfection.

When cleaning, ensure that all residues are removed. Never use any substances that may be poisonous when in contact with the skin or mucous membranes, or when swallowed or inhaled.

Clean the device using a damp cloth, to which you can apply a little detergent if necessary.

Do not use abrasive detergents and never immerse the unit in water.

1. Ensure that no water penetrates inside the unit!
2. Do not put the device or accessories in a dishwasher.
3. Do not touch the unit with wet hands while it is plugged in; do not allow any water to be sprayed onto the unit. The unit must be operated only when it is completely dry.
4. Do not spray any liquids into the ventilation slots. Any liquid that penetrates into the unit can damage the electrical parts or other components of the nebulizer and impair the function of the device.

## 6. DISINFECTION

Follow the steps listed below carefully in order to disinfect your nebulizer bottle and the accessories. It is advisable to disinfect the individual components after the last time they are used each day at the latest.

1. First clean the nebulizer bottle and the accessories as described under "CLEANING".
2. Use one part white vinegar with three parts distilled water. Make sure mixed solution is enough to submerge the nebulizer bottle, air tube, mouth piece and mask.
3. Wash nebulizer bottle, mouth piece and mask in warm water and a mild detergent. Then wash them in hot tap water.
4. Submerge these parts in the vinegar and water solution for thirty minutes.
5. You can also disinfect the nebulizer bottle, mouth piece, mask and air tube with a standard vaporizer. Observe the manufacturer's instructions for use for the vaporizer when doing this.

**⚠ Note:** The nebulizer bottle and mask must not be placed in boiling water!

6. Carefully dry all components using a soft cloth.
7. We recommend disinfecting the masks using a retail disinfectant.
8. Once all components are fully dry, reassemble the nebulizer bottle and place the components in a dry, sealed container.

**⚠ Note:** Ensure that all parts are thoroughly dried during, and in particular after cleaning, otherwise there is an increased risk of bacterial growth.

### Drying

Place the individual components on a dry, clean and absorbent surface and allow to dry completely (for at least 4 hours).

## 7. REPLACEMENT PARTS & OPTIONAL ACCESSORIES

Replacement parts are available through your local retailer.

Optional accessories include child mask, nose piece, baby pacifier inhaler. The pacifier inhaler allows baby to receive aerosol treatment easily through normal breath. (Figure 8)

The air-flow control device can control the release of air-flow, so the medicine can be saved when the treatment pauses (Figure 8).


Fig. 8

### Replacement of the air filter

Check the filter after every 10<sup>th</sup> use. If the filter is very contaminated or blocked, replace the used filter. If the filter has become wet, it must also be replaced with a new one.

1. Remove the filter cap from the unit (Figure 9).
2. Replace the used filter with a new one (Figure 10).
3. Re-insert the filter cap in its position (Figure 11).

**⚠ Note:** Do not try to clean and re-use the used filter.

Only use original filters supplied by the manufacturer. Otherwise, your nebulizer may become damaged or effective therapy can no longer be guaranteed.

Never operate the device without a filter.


Fig. 9


Fig. 10


Fig. 11

## 8. TROUBLESHOOTING

Problem/Question	Possible Cause/Remedy
The nebulizer bottle produces little or no aerosol.	<p>1. Too much or too little medication in the nebulizer bottle. Minimum: 2 ml, maximum: 8 ml, recommended: 5 ml</p> <p>2. Check that the nozzles are not blocked. If necessary, clean the nozzle (e.g. by rinsing). Use the nebulizer bottle again. <b>Caution:</b> Only poke implements through the fine holes from the under-side of the nozzle.</p> <p>3. The nebulizer bottle is not held in an upright position.</p> <p>4. The medicine solution is unsuitable for nebulizing (e.g. too thick). <b>The medicine solution should be specified by the doctor.</b></p>
Output is too low.	Air tube kinked, filter blocked, too much inhalation solution.
Which medicines are suitable for inhalation?	<b>Consult your doctor.</b> Essentially all medicines that are suitable and approved for device-based inhalation can be inhaled.
Some inhalation solution remains in the nebulizer bottle.	This is normal and occurs for technical reasons. Stop inhalation if you hear a notable difference in the noise of the nebulizer bottle.
What special steps should be taken for babies and children?	<p>1. In babies, the mask should cover the mouth and nose to guarantee effective inhalation.</p> <p>2. In children, the mask should also cover the mouth and nose. Nebulization next to a sleeping person is not suitable because insufficient medicine reaches the lungs. <b>Note:</b> Children should only use the device with help and under supervision of an adult. Never leave a child alone with the nebulizer.</p>
Inhalation with the mask takes longer.	This is for technical reasons. Less medicine is inhaled per breath through the mask holes than using the mouthpiece. The aerosol is mixed with ambient air through the holes in the mask.
Why do I need to replace the nebulizer at regular intervals?	<p>There are two reasons for this:</p> <p>1. To guarantee a therapeutically effective particle spectrum, the nozzle holes should not exceed a certain diameter. Due to the thermal and mechanical stresses, the plastic is subjected to a certain amount of wear. The nebulizer bottle is particularly sensitive. This can also cause changes to the droplet composition of the aerosol, which has a direct effect on the effectiveness of the treatment.</p> <p>2. You are also recommended to replace the nebulizer bottle on a regular basis for hygiene reasons.</p>
Does each person need their own nebulizer bottle?	This is absolutely necessary for hygiene reasons.

## 9. SPECIFICATIONS

Power	AC 230V 50 Hz
Power Consumption	180 VA
Medication Capacity	5ml
Particle Size	0.5 to 10 µm
MMAD	5 µm
Sound Level	around 60 dB
Average Nebulization Rate	0.2 ml/min.
Compressor Pressure Range	25 to 42 psi (0.17 to 0.29 Mpa)
Operating Pressure Range	8 to 13 psi (0.056 to 0.09 Mpa)
Liter Flow Range	5.0 ~ 8.0 Lpm
Operating Temperature Range	10 to 40 °C (50 to 104 °F)
Operating Humidity Range	10 to 95% RH
Storage Temperature Range	-21 to 70 °C (-13 to 158 °F)
Storage Humidity Range	10 to 95 % RH
Dimension (L x W x H)	200 × 150 × 92 mm
Weight	1.324 kg
Standard Accessories	Nebulizer, Air Tube, Mouth piece, adult mask, 5 filters. It can be equipped with a children's mask
Optional Accessories	Baby pacifier inhalator w/connector, child mask, medicine saving device

Consider instructions for use

Application part of Type BF

Attention!

Consider instructions for use

Disposal for separate collection

Vega Technologies Inc.  
11F-13, 100 Chang Chun Road,  
Taipei, Taiwan

Biolife, Ltd.  
203, Mesogion Av. P.C. 115 25 Athens,  
Greece.(VAT NO: EL 999725460)

ISO9001 & ISO13485 Certified  
This device complies with EU EMC standard EN60601-1-2

## 10. EMC INFORMATION

### Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Emissions	Compliance	Electromagnetic environment—guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RP emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

### Recommended separation distances between portable and mobile RF communications equipment and the IH25/1

The IH25/1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the IH25/1 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IH25/1 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3,8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

#### Notes:

- At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 KV contact ±8 KV air	±6 KV contact ±8 KV air	Floors should be wood, concrete or tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 KV for power supply lines	±2 KV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 KV line(s) and neutral	±1 KV line(s) and neutral	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles >5 % UT (>95 % dip in UT) for 5s	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles >5 % UT (>95 % dip in UT) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Mains power quality should be that of a typical commercial or hospital environment.

**Note:** UT is the a.c. mains voltage prior to application of the test level.

### Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range.** Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

#### Notes:

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

\*. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can not be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the IH25/1 is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

\*\* Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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